


2025/9/26



Robot 2 position for double Cd.

Zm: 292012


Zout: 199560

Nb.	Position	Block path	Open path.
1	-27000	1	2-3 > IFG PS2
2	0	2	1-3 > IFG PS1 and PS2
3	94K	3	1-2 > IFG PS2
4	48500	1-2	3 > <del>IFG PS</del>
5	72500	2-3	1 >
6	25700	3-1.	2 >

Run IFG-s for each double Cd blocker position

iPg2-3p20s-B2-26Sep1859-inf


iPg1-Bp20s-B3-27Sep1536-inf



27/09/2025

Run camera measurements of the new aperture

5mm



7.5mm

$\gamma = 18.5$

$z = 15$

iPg2-3p200s-B1-27Sep2010-inf

iPg1-3p200s-B3-28Sep1327-inf

2


28/09/2025

Scan different aperture positions and size

3 mm



5 mm



Aperture\_scan - ifg1 - 3p30s - 28 Sep 1825.inf  
Should be  
2

Repeat scan during night with smaller aperture

3 mm



3 mm

Aperture\_scan - ifg1 - B - 3p40s - 28 Sep inf

Aperture\_scan - ifg2 - B - 3p40s - 29 Sep 0307.inf

29/09/2025

Aperture\_scan - ifg1 - B3 - 3p40s - 29 Sep 1122.inf

Try aperture pos

3 mm



3 mm

 $z = 10$  $y = 18$ 

- ①

3 mm




3 mm

 $z = 13.5$  $y = 18$ 

- ②

①, ② have the same intensities tendency.



when consider  $I_{II}$  is 100%

$I_{II}$  is  $\approx 130\%$


$I_I$  is  $\approx 110\%$

robot\_diag\_double - Cd - 90deg - fat - 29 Sep 1846.inf



change phase shifter with 5mm quarts

3




previous phase shifter has 55mm high  
(quartz or sapphire)

→ change the height of phase shifter holder -


roughly 1° gives 1 period

double-Cd absorber setting back




No	Block Path	position
1	I	-21,000
2	II	5,000
3	II	95,000

Aperture Size : 3x3mm



Z




position sweep scan

6 position

x

→ 24 and 25

4 IFB




4 30/09/2025


~~ifg~~ Aperture\_scan\_ifg2\_newPS\_B1-2p50s-29Sep2253

Aperture\_scan\_ifg1\_newPS\_B3-2p50s-30Sep0524.inf

APERTURE SETTINGS:



$y = 19$      $z = 8.5$



~73%    ~64%    ~64%

IFGs look unstable despite good contrast, changed  
~~room~~ ~~temperature~~ temperature to 23°C  
Beam aperture scans during night.

Aperture\_scan\_ifg2\_newPS\_B1-2p30s-30Sep1837.inf

Aperture\_scan\_ifg1\_newPS\_B2-2p30s-01Oct0836.inf


01/11/202

IFGs are still unstable, suspected air flow in  
IFG box.

OPEN BOX LID → ifg1\_newPS-6p10s-B2-01Oct1056.inf (BAD)

CLOSED LID → ifg1\_newPS-6p10s-B2-01Oct1133.inf (GOOD)


The motor moving the PSs is a source of heat, moved  
outside the box.



(TOOK PICTURE)


5

Stability improved, running aperture scan during the night



Aperture-scan-fig 2-new PS-B1-6p10s-01 Oct 1928.inf

Aperture-scan-fig 1-new PS-B3-6p10s-01 Oct 2236.inf




$z = 11$   $y = 16.5$


02/10/2025

Scan paths intensities

robot-disg-double-Cl-test-02 Oct 1000.inf



Insert Indium to adjust intensities




free space for robot in the middle!

(INFO INDIUM: CYCLE 197 3-16-19 PAGE 175)

6

robot\_diag\_double Cd-90deg-test-02 Oct 1000.inf



ifg2-newPS-6p20s-B1-02 Oct 1244.inf

ifg1-newPS-6p20s-B2-02 Oct 1317.inf


ifg2-newPS-6p20s-B2-02 Oct 1348.inf

ifg1-newPS-6p20s-B3-02 Oct 1421.inf

ifg1-newPS-6p20s-B3-03 Oct 0949.inf

03/10/2025

Select PS angles ranges


PS1  $\rightarrow [-1, 1]$ PS2  $\rightarrow [-3, -1]$ 

ifg2-newPS-2p30s-B1-03 Oct 1202.inf C ~ 70%

ifg1-newPS-2p30s-B2-03 Oct 1228.inf C ~ 71%

ifg2-newPS-2p30s-B2-03 Oct 1252.inf C ~ 72%

ifg1-newPS-2p30s-B3-03 Oct 1318.inf C ~ 65%




ifg2-newPS-2p30s-B1-03 Oct 1357.inf C ~ 70%

ifg1-newPS-2p30s-B2-03 Oct 1429.inf C ~ 71%

ifg2-newPS-2p30s-B2-03 Oct 1454.inf C ~ 71%

ifg1-newPS-2p30s-B3-03 Oct 1520.inf C ~ 68%






ifg2 - new PS - 2p30s - B1 - 03 Oct 1611.inf C ~ 69%

ifg1 - new PS - 2p30s - B2 - 03 Oct 1<sup>637</sup>~~600~~.inf C ~ 69%

ifg2 - new PS - 2p30s - B2 - 03 Oct 1702.inf C ~ 70%

ifg1 - new PS - 2p30s - B3 - 03 Oct 1728.inf C ~ 65%



Aperture - scan - ifg<sup>2</sup>~~1~~ - new<sup>PS</sup>~~PS~~ - B1 - 2p30s - 03 Sep 1911.inf

Aperture - scan - ifg1 - new PS - B3 - 2p30s - 04 Oct 1232.inf

04/10/2025

Aperture - scan - ifg2 - new PS - B1 - 2p30s - 04 Oct 1906.inf

~~Aperture~~ - scan - ifg1 - new PS - B3 - 2p30s - 04 Oct 2301.inf

05/10/2025

z = 10 y = 16.5

ifg2 - new PS - 2p30s - B1 - 05 Oct 1433.inf C ~ 65%

ifg1 - new PS - 2p30s - B2 - 05 Oct 1500.inf C ~ 72%

ifg2 - new PS - 2p30s - B2 - 05 Oct 1527.inf C ~ 72%

ifg1 - new PS - 2p30s - B3 - 05 Oct 1555.inf C ~ 66%

z = 11.5 y = 16.5

ifg2 - new PS - 2p30s - B1 - 05 Oct 1624.inf C ~ 68%

ifg1 - new PS - 2p30s - B2 - 05 Oct 1651.inf C ~ 68%

ifg2 - new PS - 2p30s - B2 - 05 Oct 1718.inf C ~ 68%

ifg1 - new PS - 2p30s - B3 - 05 Oct 1745.inf C ~ 63%

8

Z = 10.5 Y = 16.5

ifg newPS-2p30s-B1-05Oct1814.inf C ~ 65%

Z = 11 Y = 16.5

no intensity after ifg2-newPS-2p30s-B1-05Oct2104 → ok  
 06/10/2025 rocking-05Oct2121 → no intensity

Robot bumped into IFM box cover and changed Col stimulation. Col bumped into IFM, no damage from measurements with camera. Changed ZOUT-DOUBLE Col → 207413

ifg1-newPS-2p120s-B3-06Oct1208.inf/tiF

ifg2-newPS-3p120s-B2-06Oct1340.inf/tiF

Updated Robot diagonal positions:

BLOCK	POSITION (DIAG.)	NEW LABELS (SEE PAGE 10)
1	-23000	3
2	0	2
1-3	26000	1-3
1-2	46000	3-2
2-3	70000	1-2
3	90000	1

Measured single path intensities after placing Indium

PATH COUNTS (60s)

1 1547

2 1545

3 1525

Aperture-sc3h-ifg2-newPS-B1-6p20s-06Oct2117.inf




7.10.25 8:30 reactor power 47 MW  $\rightarrow$  55 MW

$\delta$	FWHM scan	in $z_s$
-1.2	0,000 642	953
-1.15	530	1108
-1.10	477	1161
-1.05	59	996
-1.12	481	1160

Aperture - scan - ifg 1 - new PS - B2 - 6p10s - 07 Oct 1119.inf

Aperture - scan - ifg 1 - new PS - B3 - 6p10s - 07 Oct 1532.inf



Select aperture position and PSs range:

$Y=15$   $Z=11$  PS1  $[-1.5, 0.5]$  PS2  $[-1, 1]$

Test of measurement of weak value IFGs

~~preparation~~  $w_{+1} = \frac{\langle + | \hat{\Pi}_1 | \psi \rangle}{\langle + | \psi \rangle} = \text{"WV 1"}$

$$|\psi^{x_1 x_2 x_3}\rangle = \frac{1}{\sqrt{3}} (e^{-ix_1} |1\rangle + e^{-ix_2} |2\rangle + e^{-ix_3} |3\rangle)$$

$$= \text{"psi-} x_1 - x_2 - x_3 \text{"}$$

$|1\rangle \rightarrow 0$  - BEAM

prep - psi - 0 - 0 - 0.sc  $\rightarrow |\psi^{000}\rangle = |+\rangle = \frac{1}{\sqrt{3}} (|1\rangle + |2\rangle + |3\rangle)$

ifg 1 R - new PS - 2p30s - 0 - 07 Oct 1813.inf

ifg 1 R - new PS - 2p30s - p90 - 08 Oct 1143.inf

RELATIVE STEP

PHASE SHIFT

08/10/2025


Mistake in the script "ifg1r\_newPS-2p30s.sc" fixed  
(It started with a step instead of measurement.)

Adjusted fit cosine to keep the IFBs in the same  
PS range (or almost)


## CHANGED PATHS LABELS!!!

The ~~gate~~ labeling was not consistent with previous  
experiments. Changed script names accordingly ("NL" → NEW  
LABEL)

### OLD LABELS



### NEW LABELS!



PREP

### STATE PREPARATION

SET  $\chi_2 = \chi_3$  (IFG 1 + FIT AX)

SET  $\chi_1 = \chi_2 = \chi_3$  (IFG 2 + FIT 0)

RESULT →  $|\Psi\rangle = \frac{1}{\sqrt{3}}(|11\rangle + |12\rangle + |13\rangle)$

7:52  
3.10.25 address error when moving robot diagonal

ifg1r\_newPS-2p30s-psi-1+1+1\_0\_08Oct1314.inf

ifg1r\_newPS-2p30s-psi-1-1-1-1008Oct0904.inf

Relative  
steps

$$|\Psi_{IN}\rangle = \frac{1}{\sqrt{3}}(|11\rangle - |12\rangle - |13\rangle)$$

+180  $\Delta\chi_1$



Run different measurement tests, procedure still not optimal (Script mistakes, fit failing to select right phases).

Changed measurement file structure

<sup>WEAK VALUE  $P_1$</sup>   
ifg - WV1 -  <sup>$\psi_1$</sup>  psi - +1+1+1 - ..... .inf

- ifg - WV1 ..... .dat  $\leftrightarrow \Delta x_1 = 0$
- ifg - WV1 ..... .dat  $\leftrightarrow \Delta x_1 = -\frac{\pi}{2}$
- ifg - WV1 ..... .dat  $\leftrightarrow \Delta x_1 = \frac{\pi}{2}$
- ifg - WV1 ..... .dat  $\leftrightarrow \Delta x_1 = \pi$
- ifg - WV1 ..... .dat  $\leftrightarrow$  path intensity (2 $\rightarrow$ 1 $\rightarrow$ 3)

ifg - WV1 - psi - +1+1+1 - 09 Oct 1515 .inf } Wrong PS pos.  
ifg - WV1 - psi - +1+1+1 - 08 Oct 1951 .inf }

10/10/25



Measurement tests and adjustments, some are good

13/10/25



$\chi^2 = 16, z = 11$

Contrast - loops - B1 - 2 - 3 - 13 Oct 1020 .inf



B1

$C \approx 0.62$



B2

$C \approx 0.74$



B3

$C \approx 0.69$

20:36 Address error on Rohd Diagonal




12

13/10  $\rightarrow$  16/10 Measurements, optimization and attempts to prepare  $|\psi\rangle = \frac{1}{\sqrt{3}} (|1\rangle + e^{-i\frac{\pi}{3}}|2\rangle + e^{i\frac{2\pi}{3}}|3\rangle)$

Changed PSs movement direction of  $x_2$  scans

$x_1$  SCAN



$x_2$  SCAN



$x_3$  SCAN



16/10  $\rightarrow$  20/10

Performed different measurements, LABVIEW often crashes or fits do not converge. List of good and bad measurements present in the python scripts. We decided to change measurement method to ~~avoid~~ avoid fitting

**NEW METHOD:** EXTRACT PERIOD FROM FITS AND TAKE EXACTLY 16 POINTS PER PERIOD. DATA IS OBTAINED BY SHIFTING THE POINTS.

finish 25.10. 7:00